

Rat Cyclic AMP-responsive element-binding protein 1(CREB1) ELISA Kit



Catalog No: #EK20889

Orders: order@signalwayantibody.com

Package Size: #EK20889-1 48T #EK20889-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cyclic AMP-responsive element-binding protein 1(CREB1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:31.2-2000pg/mL
Sensitivity:11.6pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only