

Mouse Growth arrest and DNA damage-inducible protein GADD34(GADD34) ELISA Kit

Catalog No: #EK20879

Orders: order@signalwayantibody.com

Package Size: #EK20879-1 48T #EK20879-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Growth arrest and DNA damage-inducible protein GADD34(GADD34) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL
Sensitivity:11.6pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only