

Mouse Jun dimerization protein 2(JDP2) ELISA Kit

Catalog No: #EK20895



Package Size: #EK20895-1 48T #EK20895-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Jun dimerization protein 2(JDP2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.058ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only