

## Collagen II Mouse mAb

Catalog No: #62081

Package Size: #62081-1 50ul #62081-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Collagen II Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous Collagen II protein
Immunogen Description	Synthetic Peptide of Collagen II
Target Name	COL2A1
Other Names	Collagen alpha-1 ; II ; chain ; Alpha-1 type II collagen ; [Cleaved into: Collagen alpha-1 ; II ; chain ; Chondrocalcin
Uniprot	P02458
GeneID	1280
Calculated MW	142kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

IHC 1:50-300; IF 1:50-200

## Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

---

Note: This product is for in vitro research use only