

eIF4G Rabbit pAb

Catalog No: #62083

Package Size: #62083-1 50ul #62083-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	eIF4G Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	eIF4G Polyclonal Antibody detects endogenous levels of eIF4G protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human eIF4G.
Target Name	EIF4G1
Other Names	EIF4G1 ; EIF4F ; EIF4G ; EIF4GI ; Eukaryotic translation initiation factor 4 gamma 1 ; eIF-4-gamma 1 ; eIF-4G 1 ; eIF-4G1 ; p220
Uniprot	Q04637
GeneID	1981
Calculated MW	175kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000; IHC 1:100-1:300; IF 1:200-1:1000; ELISA 1:5000

Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using eIF4G Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of HeLa cells, using eIF4G Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only