

## HECTD3 Antibody

Catalog No: #62200

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | HECTD3 Antibody   |
| Brief Description     | Rabbit Polyclonal   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB  |
| Species Reactivity    | Human; Mouse  |
| Specificity           | The antibody detects endogenous level of total HECTD3 protein   |
| Immunogen Type        | Peptide   |
| Immunogen Description | Synthesized peptide derived from human HECTD3   |
| Target Name           | HECTD3;   |
| Uniprot               | Q5T447;   |
| Calculated MW         | Predicted band size: 95 kDa   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                       |

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only