

## VATL Antibody

Catalog No: #62222

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VATL Antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB,IHC
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous level of total VATL protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human N-terminal VATL
Conjugates	Unconjugated.
Target Name	ATP6V0C
Uniprot	P27449
GeneID	527
Calculated MW	17kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000; IHC 1:50-200

Note: This product is for in vitro research use only