

rat IgA Antibody FITC Conjugated

Catalog No: #C00792F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	rat IgA Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Rt
Immunogen Description	Full length plasma protein from rat
Conjugates	FITC
Target Name	rat IgA
Other Names	A2m marker; FLJ14473; FLJ35065; FLJ35500; FLJ36402; FLJ39698; FLJ40001; FLJ41548; FLJ41552; FLJ41789; FLJ43248; FLJ43594; FLJ44293; FLJ46028; FLJ46621; FLJ46724; FLJ46811; FLJ46824; Ig alpha 1 chain C region; Ig alpha 2 chain C region; IgA1; Igh2; IGHA 1; IGHA 2; IGHA1; IGHA2; Immunoglobulin alpha 1
Excitation Emission	494nm 518nm
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

rat IgA

Note: This product is for in vitro research use only