

HIV1 gp41/gp160 antibody

Catalog No: #62155

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HIV1 gp41/gp160 antibody
Brief Description	Rabbit Polyclonal
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB, IHC, IF, ICC/IF, ELISA
Species Reactivity	HIV2
Immunogen Description	KLH conjugated synthetic peptide derived from HIV2 gp41 + gp160: 621-700/858
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000, IHC-P 1:100-500, IHC-F 1:100-500, IF 1:100-500, ICC/IF 1:100-500, ELISA 1:5000-10000

Note: This product is for in vitro research use only