

Dopamine Receptor D1 Antibody PE Conjugated

Catalog No: #C00816P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dopamine Receptor D1 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt, ICC, IF
Species Reactivity	Hu Rt
Immunogen Description	KLH conjugated synthetic peptide aa 2-23 446 derived from human Dopamine Receptor D1
Conjugates	PE
Target Name	Dopamine Receptor D1
Other Names	DADR; DRD1A; D(1A) dopamine receptor; Dopamine D1 receptor; DRD1
Accession No.	Swiss-Prot#P21728NCBI Gene ID1812
Uniprot	P21728
GeneID	1812;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Dopamine receptor whose activity is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only