

## CK8 Antibody FITC Conjugated

Catalog No: #C00927F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CK8 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 483 derived from human CK8
Conjugates	FITC
Target Name	CK8
Other Names	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2; Keratin, type II cytoskeletal 8; Cytokeratin-8; Keratin-8; Type-II keratin Kb8; KRT8
Accession No.	Swiss-Prot#P05787NCBI Gene ID3856
Uniprot	P05787
GeneID	3856;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test IF=1:50-200

## Background

Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Note: This product is for in vitro research use only