

Myelin PLP Antibody FITC Conjugated

Catalog No: #C00944F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Myelin PLP Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | ICC,IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 20-70 277 derived from human Myelin proteolipid protein |
| Conjugates | FITC |
| Target Name | Myelin PLP |
| Other Names | PLP; PMD; HLD1; MMPL; SPG2; GPM6C; PLP DM2; Myelin proteolipid protein; Lipophilin; PLP1 |
| Accession No. | Swiss-Prot#P60201NCBI Gene ID5354 |
| Uniprot | P60201 |
| GeneID | 5354; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Extracellular, Cell membrane |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Background

This is the major myelin protein from the central nervous system. It plays an important role in the formation or maintenance of the multilamellar structure of myelin.

Note: This product is for in vitro research use only