

NLGN4Y Antibody FITC Conjugated

Catalog No: #C00950F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NLGN4Y Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human NLGN4Y
Conjugates	FITC
Target Name	NLGN4Y
Other Names	Y-linked; Neuroligin Y; Neuroligin-4; NLGN4Y; NLGNY_HUMAN.
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

NLGN4Y

Note: This product is for in vitro research use only