

## EEA1 Antibody FITC Conjugated

Catalog No: #C01053F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EEA1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 1280-1325 1411 derived from human EEA1
Conjugates	FITC
Target Name	EEA1
Other Names	MST105; ZFYVE2; MSTP105; Early endosome antigen 1; Endosome-associated protein p162; Zinc finger FYVE domain-containing protein 2; EEA1
Accession No.	Swiss-Prot#Q15075NCBI Gene ID8411
Uniprot	Q15075
GeneID	8411;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

Binds phospholipid vesicles containing phosphatidylinositol 3-phosphate and participates in endosomal trafficking.

Note: This product is for in vitro research use only