

HB9 HLXB9 Antibody PE Conjugated

Catalog No: #C01115P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HB9 HLXB9 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 250-300 401 derived from human HLXB9
Conjugates	PE
Target Name	HB9 HLXB9
Other Names	HB9; HLXB9; SCRA1; HOXHB9; Motor neuron and pancreas homeobox protein 1; Homeobox protein HB9; MNX1
Accession No.	Swiss-Prot#P50219NCBI Gene ID3110
Uniprot	P50219
GeneID	3110;
Excitation Emission	480,565nm 578nm
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Putative transcription factor involved in pancreas development and function.

Note: This product is for in vitro research use only