

NG Neurogranin Antibody HRP Conjugated

Catalog No: #C01205H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NG Neurogranin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F,ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-78 78 derived from human NG Neurogranin
Conjugates	HRP
Target Name	NG Neurogranin
Other Names	RC3; hng; Neurogranin; Ng; NRGN
Accession No.	Swiss-Prot#Q92686NCBI Gene ID4900
Uniprot	Hu Ms Rt
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Acts as a third messenger" substrate of protein kinase C-mediated molecular cascades during synaptic development and remodeling. Binds to calmodulin in the absence of calcium (By

similarity).,Q92686,4900,;https://www.uniprot.org/uniprot/Q92686,https://www.ncbi.nlm.nih.gov/gene/?term=4900

33977,33977,C01207F,,TRIP13 16E1BP Antibody FITC Conjugated,Rabbit,Polyclonal,,IgG,Purified by Protein A.,ICC,IF

References

ICC=1:50-200,IF=1:50-200,

Note: This product is for in vitro research use only