

SEMA4F Antibody FITC Conjugated

Catalog No: #C01238F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SEMA4F Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 770 derived from human SEMA4F
Conjugates	FITC
Target Name	SEMA4F
Other Names	S4F; SEMAM; SEMAW; M-SEMA; PRO2353; m-Sema-M; Semaphorin-4F; Semaphorin-M; Sema M; Semaphorin-W; Sema W; SEMA4F
Accession No.	Swiss-Prot#O95754NCBI Gene ID10505
Uniprot	O95754
GeneID	10505;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Has growth cone collapse activity against retinal ganglion-cell axons.

Note: This product is for in vitro research use only