

SEMA4F Antibody HRP Conjugated

Catalog No: #C01238H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SEMA4F Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F,ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 770 derived from human SEMA4F
Conjugates	HRP
Target Name	SEMA4F
Other Names	S4F; SEMAM; SEMAW; M-SEMA; PRO2353; m-Sema-M; Semaphorin-4F; Semaphorin-M; Sema M; Semaphorin-W; Sema W; SEMA4F
Accession No.	Swiss-Prot#O95754NCBI Gene ID10505
Uniprot	O95754
GeneID	10505;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Has growth cone collapse activity against retinal ganglion-cell axons.

Note: This product is for in vitro research use only