

Bid Antibody FITC Conjugated

Catalog No: #C01290F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bid Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 142-192 195 derived from mouse Bid
Conjugates	FITC
Target Name	Bid
Other Names	AI875481; AU022477; 2700049M22Rik; BH3-interacting domain death agonist; p22 BID; BID
Accession No.	Swiss-Prot#P70444NCBI Gene ID12122
Uniprot	P70444
GeneID	12122;
Excitation Emission	494nm 518nm
Cell Localization	Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Induces caspases and apoptosis. Counters the protective effect of Bcl-2. The major proteolytic product p15 BID allows the release of cytochrome c.

Note: This product is for in vitro research use only