

AGEs Antibody HRP Conjugated

Catalog No: #C01328H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AGEs Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC
Species Reactivity	Hu Ms Rt
Immunogen Description	AGEs
Conjugates	HRP
Target Name	AGEs
Other Names	Serum albumin; BSA; ALB
Accession No.	Swiss-Prot#P02769NCBI Gene ID280717
Uniprot	P02769
GeneID	280717;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 8% of all plasma zinc.

Note: This product is for in vitro research use only