

TFEB Antibody


Catalog No: #25476

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TFEB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	TFEB antibody was raised against an 18 amino acid peptide near the amino terminus of human TFEB.
Target Name	TFEB
Other Names	Transcription factor EB, TCFEB, BHLHE35, ALPHATFEB
Accession No.	Swiss-Prot:P19484Gene ID:7942
Uniprot	P19484
GeneID	7942;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of TFEB in A549 cell lysate with TFEB antibody at (A) 1 and (B) 2 ug/mL.

Note: This product is for in vitro research use only