

## PRC1 Antibody

Catalog No: #25481


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PRC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Specificity	At least three alternatively spliced transcript variants encoding distinct isoforms have been observed.
Immunogen Type	Peptide
Immunogen Description	19 amino acid peptide near the center of human PRC1.
Target Name	PRC1
Other Names	Protein Regulator of Cytokinesis 1, ASE1, MGC1671, MGC3669
Accession No.	Swiss-Prot:O43663Gene ID:9055
Uniprot	O43663
GeneID	9055;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Images



Western blot analysis of PRC1 in human skeletal muscle tissue lysate with Prc1 antibody at (A) 0.5 ug/ml and (B) 1 ug/ml.

Note: This product is for in vitro research use only