

## CLEC7A Antibody


Catalog No: #25485

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CLEC7A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Specificity	Multiple isoforms of CLEC7A are known to exist.
Immunogen Type	Peptide
Immunogen Description	17 amino acid peptide near the center of human CLEC7A.
Target Name	CLEC7A
Other Names	C-type lectin domain family member 7 isoform A, BGR, DECTIN1, CLECSF12, CANDF4
Accession No.	Swiss-Prot:Q9BXN2Gene ID:64581
Uniprot	Q9BXN2
GeneID	64581;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Images



Western blot analysis of CLEC7A in rat spleen tissue lysate with CLEC7A antibody at (A) 1 and (B) 2 ug/.

Note: This product is for in vitro research use only