

ABCA7 Antibody


Catalog No: #25499

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABCA7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB IHC IF
Species Reactivity	Hu
Specificity	ABCA7 antibody is human specific. At least two alternatively spliced transcript variants encoding distinct isoforms have been observed.
Immunogen Type	Peptide
Immunogen Description	ABCA7 antibody was raised against an 18 amino acid peptide near the amino terminus of human ABCA7.
Target Name	ABCA7
Other Names	ATP-binding Cassette Sub-family A (ABC1) member 7, ABC1, ABC Transporter 1, ATP-binding Cassette 1, ATP-binding Cassette Transporter 1, CERP
Accession No.	Swiss-Prot:Q8IZY2Gene ID:10347
Uniprot	Q8IZY2
GeneID	10347;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.


Images



Western blot analysis of ABCA7 in 293 cell lysate with ABCA7 antibody at (A) 1 and (B) 2 ug/ml.



Immunofluorescence of ABCA7 in 293 cells with ABCA7 antibody at 20 $\mu\text{g/mL}$.



Immunohistochemistry of ABCA7 in human spleen tissue with ABCA7 antibody at 5 $\mu\text{g/mL}$.

Note: This product is for in vitro research use only