

BCL6 Antibody

Catalog No: #25523

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BCL6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Specificity	Three alternatively spliced transcript variants have been observed
Immunogen Type	Peptide
Immunogen Description	BCL6 antibody was raised against an 18 amino acid peptide near the center of human BCL6.
Target Name	BCL6
Other Names	B-cell lymphoma 5 protein, B-cell lymphoma 6 protein, BCL5, BCL-6, BCL6A, LAZ3, Zinc finger and BTB domain-containing protein 27, ZBTB27, Zinc finger protein 51, ZNF51
Accession No.	Swiss-Prot:P41182Gene ID:604
Uniprot	P41182
GeneID	604;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 78kd

Images



Western blot analysis of BCL6 in rat lung tissue lysate with Bcl6 antibody at 1 ug/ml in (A) the absence and (B) the presence of blocking peptide

Note: This product is for in vitro research use only