

## SCUBE3 Antibody

Catalog No: #25528

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SCUBE3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms
Specificity	At least two isoforms of SCUBE3 are known to exist. SCUBE3 antibody is predicted to not cross-react with other SCUBE members
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the center of human SCUBE3.
Target Name	SCUBE3
Other Names	Signal peptide CUB domain EGF-like 3, CEGF3
Accession No.	Swiss-Prot:Q8IX30Gene ID:222663
Uniprot	Q8IX30
GeneID	222663;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Application Details

Predicted MW: 109kd

## Images



Western blot analysis of SCUBE3 in mouse kidney tissue lysate with SCUBE3 antibody at (A) 1 and (B) 2 ug/ml.

Note: This product is for in vitro research use only