

GABRR1 Antibody HRP Conjugated

Catalog No: #C01718H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GABRR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide aa 210-260 479 derived from human GABRR1
Conjugates	HRP
Target Name	GABRR1
Other Names	Gamma-aminobutyric acid receptor subunit rho-1; GABA(A) receptor subunit rho-1; GABA(C) receptor; GABRR1
Accession No.	Swiss-Prot#P24046NCBI Gene ID2569
Uniprot	P24046
GeneID	2569;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA benzodiazepine receptor and opening an integral chloride channel. Rho-1 GABA receptor could play a role in retinal neurotransmission.

Note: This product is for in vitro research use only