

JARID1C Monoclonal Antibody

Catalog No: #27002

Package Size: #27002-1 50ul #27002-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	JARID1C Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms Rt D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment(Human)
Target Name	JARID1C
Other Names	lysine-specific demethylase 5C, <input type="checkbox"/> protein SmcX, <input type="checkbox"/> selected cDNA on X, <input type="checkbox"/> Smcy homolog, X-linked, <input type="checkbox"/> Smcx homolog, X chromosome, <input type="checkbox"/> histone demethylase JARID1C, <input type="checkbox"/> JmjC domain-containing protein SMCX, <input type="checkbox"/> Jumonji/ARID domain-containing protein 1C
Accession No.	UniProt:P41229 NCBI Gene ID:8242
Uniprot	P41229
GeneID	8242;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 220kd


Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected JARID1C on 40ug of cell lysates

1:150 dilution from a previous lot detected JARID1C in Hela



Note: This product is for in vitro research use only