

## UHRF1(N-term) Monoclonal Antibody

Catalog No: #27003

Package Size: #27003-1 50ul #27003-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	UHRF1(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 1aa-408aa of human UHRF1(N-terminus).
Target Name	UHRF1
Other Names	E3 ubiquitin-protein ligase UHRF1, □huNp95, □nuclear protein 95, □RING finger protein 106, □transcription factor ICBP90, □nuclear zinc finger protein Np95, □inverted CCAAT box-binding protein of 90 kDa
Accession No.	UniProt:Q96T88 NCBI Gene ID:29128
Uniprot	Q96T88
GeneID	29128;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details


Predicted MW: 90kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected UHRF1 on 40ug of cell lysates



HeLa cells using anti-UHRF1 (N-terminus)antibody diluted  
1:200

---

Note: This product is for in vitro research use only