

## TFIIB Monoclonal Antibody

Catalog No: #27007

Package Size: #27007-1 50ul #27007-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	TFIIB Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC IHC
Species Reactivity	Hu Ms P Rt D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant full length protein (316aa) of human GTF2B.
Target Name	TFIIB
Other Names	transcription initiation factor IIB <input type="checkbox"/> S300-II <input type="checkbox"/> general transcription factor TFIIB <input type="checkbox"/> RNA polymerase II transcription factor IIB
Accession No.	UniProt:Q00403 NCBI Gene ID:2959
Uniprot	Q00403
GeneID	2959;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

## Application Details


Predicted MW: 37kd

Western blotting: 1:1000


## Images



1:1000 dilution of this antibody detected TBL1Y on 40ug of cell lysates



1:150 dilution from a previous lot detected TBL1Y in Hela



Paraffin-embedded human colon using anti- GTF2B diluted 1:50-1:100

---

Note: This product is for in vitro research use only