

SOX2 Monoclonal Antibody

Catalog No: #27015

Package Size: #27015-1 50ul #27015-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SOX2 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu P Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant full length protein (319aa) of mouse Sox2.
Target Name	SOX2
Other Names	transcription factor SOX-2 □ transcription factor SOX2 □ SRY-related HMG-box gene 2
Accession No.	UniProt:P48431 NCBI Gene ID:6657
Uniprot	P48431
GeneID	6657;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 42kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected Sox2 on 40ug of cell lysates

Note: This product is for in vitro research use only