

SOX2 Monoclonal Antibody

Catalog No: #27015

Package Size: #27015-1 50ul #27015-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | SOX2 Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2b |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu P Ms Rt |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant full length protein (319aa) of mouse Sox2. |
| Target Name | SOX2 |
| Other Names | transcription factor SOX-2 □ transcription factor SOX2 □ SRY-related HMG-box gene 2 |
| Accession No. | UniProt:P48431 NCBI Gene ID:6657 |
| Uniprot | P48431 |
| GeneID | 6657; |
| Formulation | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |
| Storage | store at -20°C |

Application Details

Predicted MW: 42kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected Sox2 on 40ug of cell lysates

Note: This product is for in vitro research use only