

PRMT6 Monoclonal Antibody

Catalog No: #27017

Package Size: #27017-1 50ul #27017-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	PRMT6 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant full length protein (375aa) of human PRMT6.
Target Name	PRMT6
Other Names	protein arginine N-methyltransferase 6 □HMT1 hnRNP methyltransferase-like 6 □histone-arginine N-methyltransferase PRMT6 □heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6
Accession No.	UniProt:Q96LA8 NCBI Gene ID:55170
Uniprot	Q96LA8
GeneID	55170;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details


Predicted MW: 42kd

Western blotting: 1:1000


Images



1:1000 dilution of this antibody detected PRMT6 on 40ug of cell lysates



1:150 dilution from a previous lot detected PRMT6 in HeLa



Paraffin-embedded human breast cancer using anti- PRMT6 diluted 1: 50-1: 100

Note: This product is for in vitro research use only