

## P53 Monoclonal Antibody

Catalog No: #27021

Package Size: #27021-1 50ul #27021-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	P53 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu D Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment of human TP53 (C-terminus)
Target Name	P53
Other Names	cellular tumor antigen p53 <input type="checkbox"/> antigen NY-CO-13 <input type="checkbox"/> phosphoprotein p53 <input type="checkbox"/> p53 tumor suppressor <input type="checkbox"/> transformation-related protein 53
Accession No.	UniProt:P04637 NCBI Gene ID:7157
Uniprot	P04637
GeneID	7157;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 53kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected TP53 on 40ug of cell lysates

---

Note: This product is for in vitro research use only