

SMC1A(C-term) Monoclonal Antibody

Catalog No: #27023

Package Size: #27023-1 50ul #27023-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	SMC1A(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 2aa-351aa of human SMC1A(N-terminus).
Target Name	SMC1A
Other Names	Sb1.8
Accession No.	UniProt:Q14683 NCBI Gene ID:8243
Uniprot	Q14683
GeneID	8243;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 143kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected SMC1A on 40ug of cell lysates

Note: This product is for in vitro research use only