

## 14-3-3 Tau Antibody HRP Conjugated

Catalog No: #C01991H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	14-3-3 Tau Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human 14-3-3 Tau
Conjugates	HRP
Target Name	14-3-3 Tau
Other Names	14 3 3 protein T cell; 14 3 3 protein tau; 14 3 3 protein theta; 14 3 3 tau; 14 3 3 theta; 14-3-3 protein T-cell; 14-3-3 protein tau; 14-3-3 protein theta; 1433T_HUMAN; 1C5; HS1; KCIP1; Protein HS1; Protein kinase C inhibitor protein 1; Protein tau; tyr3 trp5 monooxygenase activation protein, theta;
Accession No.	NCBI Gene ID10971
Uniprot	P27348
GeneID	10971;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

## Background

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. Negatively regulates the kinase activity of PDPK1. This antibody is reactive with 14-3-3 Alpha Beta Gamma Delta Epsilon.

Note: This product is for in vitro research use only