

BiP/GRP78(C-term) Monoclonal Antibody

Catalog No: #27033

Package Size: #27033-1 50ul #27033-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	BiP/GRP78(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 487aa-638aa of human HSPA5(C-terminus)
Target Name	BiP/GRP78
Other Names	78 kDa glucose-regulated protein □immunoglobulin heavy chain-binding protein □endoplasmic reticulum luminal Ca(2+)-binding protein grp78
Accession No.	UniProt:P11021 NCBI Gene ID:3309
Uniprot	P11021
GeneID	3309;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 89kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected HSPA5 (C-terminus) on 40ug of cell lysates



HeLa cells using anti- HSPA5 (C-terminus) antibody diluted 1:100

Note: This product is for in vitro research use only