

USP7(C-term) Monoclonal Antibody

Catalog No: #27041

Package Size: #27041-1 50ul #27041-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	USP7(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms Rt D P Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 868aa-1093aa of human USP7 (C-terminus).
Target Name	USP7
Other Names	ubiquitin carboxyl-terminal hydrolase 7 □ubiquitin thioesterase 7 □deubiquitinating enzyme 7 □ubiquitin-specific-processing protease 7 □herpesvirus-associated ubiquitin-specific protease □Herpes virus-associated ubiquitin-specific proteas
Accession No.	UniProt:Q93009NCBI Gene ID:7874
Uniprot	Q93009
GeneID	7874;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 129kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected USP7 on 40ug of cell lysates



HeLa cells using anti- USP7 (N-terminus) antibody diluted 1:150

Note: This product is for in vitro research use only