

Albumin Monoclonal Antibody

Catalog No: #27044

Package Size: #27044-1 50ul #27044-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Albumin Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human Albumin
Target Name	Albumin
Other Names	serum albumin <input type="checkbox"/> albumin (32 AA) <input type="checkbox"/> albumin (AA 34) <input type="checkbox"/> growth-inhibiting protein 20 <input type="checkbox"/> cell growth inhibiting protein 42
Accession No.	UniProt:P02768 NCBI Gene ID:213
Uniprot	P02768
GeneID	213;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details


Predicted MW: 68kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected Albumin on 40ug of cell lysates



Paraffin--embedded human liver cancer using anti- ALBdiluted
1:50-1:100

Note: This product is for in vitro research use only