

FAS Monoclonal Antibody

Catalog No: #27054

Package Size: #27054-1 50ul #27054-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FAS Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 211aa-329aa of human CD95(C-terminus).
Target Name	FAS
Other Names	tumor necrosis factor receptor superfamily member 6 <input type="checkbox"/> Fas AMA <input type="checkbox"/> FAS 827dupA <input type="checkbox"/> CD95 antigen <input type="checkbox"/> FASLG receptor <input type="checkbox"/> apoptosis antigen 1 <input type="checkbox"/> Delta Fas/APO-1/CD95 <input type="checkbox"/> APO-1 cell surface antigen <input type="checkbox"/> apoptosis-mediating surface antigen FAS <input type="checkbox"/> tumor nec
Accession No.	UniProt:P25445 NCBI Gene ID:355
Uniprot	P25445
GeneID	355;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 45kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected FAS on 40ug of cell lysates

Note: This product is for in vitro research use only