

## Ikaros(C-term) Monoclonal Antibody

Catalog No: #27062

Package Size: #27062-1 50ul #27062-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Ikaros(C-term) Monoclonal Antibody   |
| Host Species          | Mouse  |
| Clonality             | Monoclonal   |
| Isotype               | IgG1   |
| Purification          | Purified by antigen-affinity chromatography.   |
| Applications          | WB   |
| Species Reactivity    | Hu Ms D  |
| Immunogen Type        | Recombinant protein  |
| Immunogen Description | Recombinant fragment, corresponding to amino acids 249-415aa of human IKAROS(C-terminus).  |
| Target Name           | Ikaros   |
| Other Names           | DNA-binding protein Ikaros □ CLL-associated antigen KW-6 □ Ikaros (zinc finger protein) □ ikaros family zinc finger protein 1 □ lymphoid transcription factor LyF-1 □ zinc finger protein, subfamily 1A, 1 (Ikaros |
| Accession No.         | UniProt:Q13422<br>NCBI Gene ID:10320   |
| Uniprot               | Q13422   |
| GeneID                | 10320;   |
| Formulation           | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol  |
| Storage               | store at -20°C   |

## Application Details

Predicted MW: 58kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected IKAROS on 40ug of cell lysates

---

Note: This product is for in vitro research use only