

## FOXO1 Monoclonal Antibody

Catalog No: #27063

Package Size: #27063-1 50ul #27063-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	FOXO1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 420-649aa of human FOXO1(C-terminus).
Target Name	FOXO1
Other Names	forkhead box protein O1 □ forkhead box protein O1A □ forkhead, Drosophila, homolog of, in rhabdomyosarcoma
Accession No.	UniProt:Q12778 NCBI Gene ID:2308
Uniprot	Q12778
GeneID	2308;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details


Predicted MW: 70kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected FOXO1(C-terminus) on 40ug of cell lysates



HeLa cells using anti- FOXO1(C-terminus)antibody diluted 1:150

---

Note: This product is for in vitro research use only