

AMACR(C-term) Monoclonal Antibody

Catalog No: #27068

Package Size: #27068-1 50ul #27068-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	AMACR(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 138-319aa of human AMACR (C-terminus).
Target Name	AMACR
Other Names	alpha-methylacyl-CoA racemase □ 2-methylacyl-CoA racemase
Accession No.	UniProt:Q9UHK6 NCBI Gene ID:23600
Uniprot	Q9UHK6
GeneID	23600;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 42kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected AMACR on 40ug of cell lysates



HeLa cells using anti- AMACR (C-terminus) antibody diluted 1:150

Note: This product is for in vitro research use only