

CDK5(N-term) Monoclonal Antibody

Catalog No: #27072

Package Size: #27072-1 50ul #27072-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	CDK5(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms Rt P D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 1-138aaof human CDK5(N-terminus).
Target Name	CDK5
Other Names	cyclin-dependent kinase 5 <input type="checkbox"/> TPKII catalytic subunit <input type="checkbox"/> protein kinase CDK5 splicing <input type="checkbox"/> cell division protein kinase 5 <input type="checkbox"/> serine/threonine-protein kinase PSSALRE <input type="checkbox"/> tau protein kinase II catalytic subunit
Accession No.	UniProt:Q00535 NCBI Gene ID:1020
Uniprot	Q00535
GeneID	1020;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 36kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected CDK5 on 40ug of cell lysates



HeLa cells using anti- CDK5 (N-terminus) antibody diluted 1:150

Note: This product is for in vitro research use only