

## Hexokinase-2( N-term) Monoclonal Antibody

Catalog No: #27074

Package Size: #27074-1 50ul #27074-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Hexokinase-2( N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB FC
Species Reactivity	Hu P Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 270-435aa of human HK2(N-terminus).
Target Name	Hexokinase-2
Other Names	hexokinase-2 <input type="checkbox"/> HK II <input type="checkbox"/> hexokinase type II <input type="checkbox"/> hexokinase-2, muscle <input type="checkbox"/> muscle form hexokinase <input type="checkbox"/> mi-2 autoantigen 240 kDa protein <input type="checkbox"/> zinc-finger helicase (Snf2-like)
Accession No.	UniProt:P52789 NCBI Gene ID:3099
Uniprot	P52789
GeneID	3099;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 102kd


Western blotting: 1:1000

## Images



1:3000 dilution of this antibody detected HK2 on 40ug of cell lysates

K562 cells stained with HK2 (red, 1/100 dilution), followed by FITC-conjugated goat anti-mouse IgG. Blue line histogram represents the isotype control, normal mouse IgG



---

Note: This product is for in vitro research use only