

## TAB1(N-term) Monoclonal Antibody

Catalog No: #27075

Package Size: #27075-1 50ul #27075-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


|                       |   |
|-----------------------|---|
| Product Name          | TAB1(N-term) Monoclonal Antibody  |
| Host Species          | Mouse   |
| Clonality             | Monoclonal  |
| Purification          | Purified by antigen-affinity chromatography.  |
| Applications          | WB  |
| Species Reactivity    | Hu D Ms P Rt  |
| Immunogen Type        | Recombinant protein   |
| Immunogen Description | Recombinant fragment, corresponding to amino acids of human TAB1  |
| Target Name           | TAB1  |
| Other Names           | TGF-beta-activated kinase 1 and MAP3K7-binding protein 1 □ TAK1-binding protein 1 □ transforming growth factor beta-activated kinase-binding protein 1 □ mitogen-activated protein kinase kinase kinase 7-interacting protein 1 |
| Accession No.         | UniProt:Q15750<br>NCBI Gene ID:10454  |
| Uniprot               | Q15750  |
| GeneID                | 10454;  |
| Formulation           | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol   |
| Storage               | store at -20°C  |

## Application Details

Predicted MW: 55kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected TAB1 on 40ug of cell lysates

---

Note: This product is for in vitro research use only