

## HNRNPU Monoclonal Antibody

Catalog No: #27076

Package Size: #27076-1 50ul #27076-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	HNRNPU Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Rb B P
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human HNRNPU
Target Name	HNRNPU
Other Names	Scaffold attachment factor A
Accession No.	UniProt:Q00839 NCBI Gene ID:3192
Uniprot	Q00839
GeneID	3192;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 90kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected HNRNPU on 40ug of cell lysates

---

Note: This product is for in vitro research use only