

## ARRB1 Monoclonal Antibody

Catalog No: #27077

Package Size: #27077-1 50ul #27077-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	ARRB1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human ARRB1
Target Name	ARRB1
Other Names	arrestin-beta-1
Accession No.	UniProt:P49407 NCBI Gene ID:408
Uniprot	P49407
GeneID	408;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 47kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected ARRB1 on 40ug of cell lysates

---

Note: This product is for in vitro research use only