

## Brg1 Monoclonal Antibody

Catalog No: #27090

Package Size: #27090-1 50ul #27090-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Brg1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4E11-E11-G11
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human SMARCA4 protein fragments expressed in E.coli.
Target Name	Brg1
Other Names	BRG1, SNF2, SWI2, RTPS2, BAF190, SNF2L4, SNF2LB, hSNF2b, SNF2-BETA
Accession No.	Swiss-Prot:P51532Gene ID:6597
Uniprot	P51532
GeneID	6597;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

## Application Details

Predicted MW: 235kd

Western blotting: 1:1000

## Images



Western blot detection of SMARCA4 in K562&HeLaNE cell lysates using SMARCA4 antibody (1:1000 diluted).  
 Predicted band size: 235KDa  
 Observed band size: 235KDa.

## Background

Transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. Component of the CREST-BRG1

complex, a multiprotein complex that regulates promoter activation.

---

Note: This product is for in vitro research use only