

GATA6 Monoclonal Antibody

Catalog No: #27092

Package Size: #27092-1 50ul #27092-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	GATA6 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3H3-H6-H11
Isotype	unknown
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human GATA6 protein fragments expressed in E.coli.
Target Name	GATA6
Other Names	ASD9, AVSD5, PACHD
Accession No.	Swiss-Prot:Q92908Gene ID:2627
Uniprot	Q92908
GeneID	2627;
Formulation	ascites

Application Details

Predicted MW: 60kd

Western blotting: 1:1000

Images



Western blot detection of GATA6 in HeLaNE, HeLa&293 cell lysates and using GATA6 antibody (1:200 diluted). Predicted band size: 45kDa Observed band size: 60kDa

Background

Thought to be important for regulating terminal differentiation and/or proliferation.

Note: This product is for in vitro research use only